


# Journal Pre-proof

Monitoring alkyl pyrazines in roasted hazelnuts by SHS-GC-IMS: IMS response assessment and standardization

Maria Mazzucotelli, Iuliia Khomenko, Emanuela Betta, Irene Cetto, Brian Farneti, Elena Gabetti, Andrea Cavallero, Eugenio Aprea, Franco Biasioli



PII: S0039-9140(23)00319-3

DOI: <https://doi.org/10.1016/j.talanta.2023.124568>

Reference: TAL 124568

To appear in: *Talanta*

Received Date: 22 February 2023

Revised Date: 14 April 2023

Accepted Date: 16 April 2023

Please cite this article as: M. Mazzucotelli, I. Khomenko, E. Betta, I. Cetto, B. Farneti, E. Gabetti, A. Cavallero, E. Aprea, F. Biasioli, Monitoring alkyl pyrazines in roasted hazelnuts by SHS-GC-IMS: IMS response assessment and standardization, *Talanta* (2023), doi: <https://doi.org/10.1016/j.talanta.2023.124568>.

This is a PDF file of an article that has undergone enhancements after acceptance, such as the addition of a cover page and metadata, and formatting for readability, but it is not yet the definitive version of record. This version will undergo additional copyediting, typesetting and review before it is published in its final form, but we are providing this version to give early visibility of the article. Please note that, during the production process, errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.

© 2023 Published by Elsevier B.V.

Maria Mazzucotelli – Methodology; Formal analysis; Investigation; Data Curation; Visualization; Writing - Original Draft; Writing - Review & Editing

Iuliia Khomenko – Methodology; Data Curation; Visualization; Writing - Review & Editing

Emanuela Betta – Investigation; Writing - Review & Editing

Irene Cetto – Investigation; Writing - Review & Editing

Brian Farneti - Data Curation; Visualization; Writing - Review & Editing



Elena Gabetti – Resources; Writing - Review & Editing

Andrea Cavallero – Conceptualization; Resources; Project administration; Writing - Review & Editing

Eugenio Aprea - Conceptualization; Project administration; Writing - Review & Editing

Franco Biasioli – Conceptualization; Project administration; Writing - Review & Editing


# Alkyl pyrazine determination by SHS-GC-IMS



Response assessment  
and standardization

Identification

Roasting  
monitoring



# 1 **Monitoring alkyl pyrazines in roasted hazelnuts by SHS-GC-IMS: IMS response** 2 **assessment and standardization**

3 Maria Mazzucotelli<sup>1,3</sup>, Iuliia Khomenko<sup>\*1</sup>, Emanuela Betta<sup>1</sup>, Irene Cetto<sup>1</sup>, Brian Farneti<sup>1</sup>, Elena Gabetti<sup>2</sup>,  
4 Andrea Cavallero<sup>2</sup>, Eugenio Aprea<sup>1,3</sup>, Franco Biasioli<sup>1</sup>

5 \*Corresponding Author

6 <sup>1</sup> Research and Innovation Center, Fondazione Edmund Mach, San Michele all'Adige (Trento), Italy

7 <sup>2</sup> Soremartec Italia srl, Alba (Cuneo), Italy

8 <sup>3</sup> Center for Agriculture Food Environment C3A, University of Trento, San Michele all'Adige (Trento), Italy

9

## 10 **Abstract**

11 Gas chromatography coupled with ion mobility spectrometry (IMS) is an analytical tool which is rapidly  
12 becoming widespread in the analysis of food volatiles. Despite this increasing popularity, an assessment of  
13 the IMS response for several flavor compound classes is not yet available. This study focuses on alkyl  
14 pyrazines and their determination in roasted hazelnut pastes. These Maillard reaction products are crucial to  
15 determine the aromatic profile of roasted foods and are suitable markers for industrial roasting monitoring.  
16 The instrumental response of 8 alkyl pyrazines was studied using a model matrix and a matrix matching  
17 approach. The results showed a relevant effect of the pyrazine ring substitution pattern on the concentration-  
18 response curve trends, highlighting that an external standardization of the IMS response is required to make  
19 possible relative abundance comparisons between analytes. A response standardization was therefore  
20 developed and applied to determine alkyl pyrazines in samples with different roasting intensity and  
21 geographical and botanical origin.

## 22 **Keywords**

23 SHS-GC-IMS; Alkyl pyrazines; Detector response; Standardization; Hazelnut; Roasting

## 1) Introduction

The development of Ion Mobility Spectrometry (IMS) traces back to the 1960s, and it has been used, for many applications, as a standalone analytical technique to detect compounds in trace amounts [1]. IMS is a gas phase separation technology, operating at atmospheric pressure, and based on two sequential processes i) the ionization of volatile organic compounds (VOCs) in a reaction region (the ion source region); ii) the separation of the ions in a drift region according to their mobility in an electric field through a counterflowing gas (drift gas). Being the mobility of ions inversely proportional to their collisional cross-section, the IMS was defined as an ion separator based on the size/charge [2].

Since the early stages of the IMS development, the potentialities of using it as a detector after a chromatographic separation have been explored [3]. Although the first paper reporting the coupling of IMS detection with gas chromatography separation was published in 1970 [4], only recently we are witnessing a rapid increase in the popularity of GC-IMS commercial instruments, especially for food volatile analysis, with a big interest in traceability and flavor characterization purposes [5, 6]. The hyphenation of these two techniques results in an analytical platform characterized by high sensitivity, low detection limits, and an orthogonal bidimensional separation of analytes.

The capability to detect VOCs at low concentrations is the main advantage of GC-IMS since it makes it possible to apply static headspace (SHS) sampling and to avoid preconcentration steps. This reflects in minimal sample preparation and relatively short analysis time (less than an hour). Additionally, from a hardware perspective, no vacuum system is needed, resulting in analytical platforms with reduced bench space requirements, simple maintenance, and low energy consumption. The combination among the user friendliness of the hardware, the limited sample preparation required, and the overall rapidity of analytical response makes SHS-GC-IMS a promising analytical tool to be applied both for academic research and industrial food processes monitoring [7], where the implementation of robust and simplified protocols is highly desirable.

On the other hand, the drawbacks of IMS detection must be taken into account when this technique is applied. First, the formation of multiple ionized species (the protonated monomer and the proton bound

49 dimer) from a single analyte is a peculiarity which significantly complicates the interpretability of the  
50 analytical output. As regards the peak identification, the mobility coefficient of the detected ions obtained  
51 with IMS detection is a useful information and libraries of mobility coefficients have been made available  
52 from instrument manufacturers. However, even their combination with retention indices is not sufficient for  
53 a reliable peak identification. In this case, analyte confirmation can only be accomplished by comparison with  
54 pure standards or using an additional MS detector [8]. Lastly, IMS has an overall limited dynamic range and  
55 shows linear response only for 1-2 orders of magnitude of concentration, due to the radioactive source.  
56 These hurdles might explain why only few published studies so far address the challenge of a quantitative  
57 approach [9 - 12], while most of the available literature focuses on the characterization of food matrices by  
58 SHS-GC-IMS applying untargeted fingerprinting or comparative cross-sample analysis.

59 In the vast majority of the recently published literature, the IMS is operated in positive ionization mode.  
60 Under these operative conditions, volatile compounds with proton affinity higher than water (691 kJ/mol  
61 [13]) are ionized and detected, thus providing a non-selective response to several chemical functionalities  
62 (alcohols, amines, ketones, aldehydes), similar to a flame ionization detector (FID) [14]. Most of the aroma  
63 compounds are, therefore, detectable by SHS-GC-IMS.

64 Clear explanations of the IMS response were reported for alkyl alcohols, aldehydes and ketones by Jurado-  
65 Campos *et al* [15], in particular related to the formation and stability of proton-bound clusters. However, a  
66 systematic assessment of the detector response for most of the other aroma compound classes has not yet  
67 been fully addressed, neither in terms of sensitivity, nor multiple ionized species formation. Moreover, the  
68 underestimation of differences in the IMS responses might lead to misinterpretations of the analytical  
69 output, in particular if assumptions are drawn from the relative intensities of peaks.

70

71 This study focuses on pyrazines, a class of volatile molecules with a well-established importance as aroma  
72 compounds. While pyrazines have been detected in food matrices characterized by SHS-GC-IMS analysis [16,  
73 17], there is no contribution specifically investigating the detector response to this class. They are  
74 heterocyclic compounds with a common chemical structure: a monocyclic aromatic ring with two atoms of

75 nitrogen in *para* position. The positions 2-, 3-, 5-, and 6- of the ring can be substituted with different groups.  
76 Alkyl pyrazines are found in thermally treated foods (generated *via* Maillard reaction) and fermented food  
77 (produced by microbial metabolism), while alkoxy pyrazines are present in raw vegetables as secondary plant  
78 metabolites [18, 19]. Due to their low odor-thresholds, they are crucial to determine the aromatic profile  
79 even though their concentrations in food matrices is low (ng/kg), thus providing an interesting case study for  
80 a GC-IMS targeted approach.

81 From an industrial perspective, alkyl pyrazines are particularly important because they are generated during  
82 the roasting processes, which greatly contribute to defining the aroma of roasted foods (coffee, cocoa, and  
83 roasted nuts – hazelnuts and peanuts) with characteristic aromatic notes (roasted, nutty, earthy, burnt).  
84 Heating conditions (temperature and duration) are known to influence not only the level of the formed  
85 pyrazines, but also their type [20]. It is therefore clear that the monitoring of these compounds during  
86 industrial roasting processes is an effective strategy to control the aroma development of roasted products.  
87 For these reasons, alkyl pyrazines have been proposed as roasting markers, and their ratios as roasting indices  
88 [21, 22].

89 The final aim of this study is the assessment of SHS-GC-IMS as a rapid technique to monitor alkyl pyrazines  
90 within the context of the industrial roasting of hazelnuts.

91 The research was structured in two phases: i) the evaluation of the GC-IMS response for a set of target alkyl  
92 pyrazines; ii) the development of an external calibration protocol to standardize the GC-IMS response and  
93 determine the content of alkyl pyrazines in roasted hazelnut paste samples from different geographical and  
94 botanical origins and processed with different roasting intensity.

## 95 **2) Materials and Methods**

### 96 **2.1 Chemicals and matrices**

97 2-methylpyrazine (2M) (CAS 109-08-0), 2,5-dimethylpyrazine (25DM) (CAS 123-32-0), 2,6-dimethylpyrazine  
98 (26DM) (CAS 108-50-9), 2-ethylpyrazine (2E) (CAS 13925-00-3), 2,3,5-trimethylpyrazine (TRI) (CAS 14667-55-  
99 1), 2,3,5,6-tetramethylpyrazine (TETRA) (CAS 1124-11-4), and a mixture of isomers of 2-ethyl-3(5or6)-

100 dimethylpyrazine (2E35DM – 61.4% - and 2E36DM – 38.6%) (CAS 27043-05-6) pure standards were  
101 purchased from Merck (Merck KGaA, Darmstadt, Germany). Stock solutions ( $20.00 \mu\text{g g}^{-1}$ ) were prepared  
102 solubilizing pure standards in a mix of medium chain triglycerides (MCT) (chain length ranging from 6 to 10  
103 carbons) provided by Oleon (Ertvelde, Evergem, Belgium). MCT was selected as solvent because, being  
104 nonvolatile, minimizes the interferences with the IMS detection. Furthermore, since the roasted hazelnut  
105 paste has a relevant fat content (50-60%), this oil-like solvent ensured miscibility and provided a possible  
106 suitable model matrix. All the diluted standard mixtures were prepared in the same solvent.

107 A deodorized roasted hazelnut paste (HZ deo), provided by Soremartec Italia Srl (Alba, Cuneo, Italy), was used  
108 as a matrix-matched blank sample. The removal of volatile components from the paste was obtained by a  
109 proprietary physical treatment.

## 110 **2.2 Hazelnut samples**

111 Hazelnut paste samples, supplied by Soremartec Italia Srl (Alba, Cuneo, Italy), were obtained processing raw  
112 kernels (*Corylus avellana* L.) from different geographical and botanical origins: Tonda Gentile Romana  
113 monocultivar hazelnuts (RO) from Lazio region (Italy), Tonda Gentile delle Langhe monocultivar hazelnuts  
114 (TGL) from Piemonte region (Italy), and Akçakoca hazelnuts (AK) from Turkey. The harvesting year was 2021  
115 for all the samples. The raw kernels were stored in sealed bags at  $4^{\circ}\text{C}$  and in controlled humidity conditions  
116 for six months prior to roasting.

117 The roasting process was carried out in a pilot scale infrared roaster at  $140^{\circ}\text{C}$ . Roasted kernels were sampled  
118 at three times (15 min, 20 min and 25 min), to determine the content of alkyl pyrazines at different roasting  
119 levels and processed to obtain the pastes. For SHS-GC-IMS analysis 1 g of paste was placed in a 20 ml HS vial.  
120 Each sample was analyzed in triplicate. All the sample vials were stored at  $-20^{\circ}\text{C}$  until the measurement.

## 121 **2.3 GC-IMS detector response curves and matrix effect assessment**

122 Concentration levels were prepared directly in the headspace vials (20 mL crimp vials), in a range between  
123  $0.1 \mu\text{g g}^{-1}$  and  $100 \mu\text{g g}^{-1}$ . Each level was prepared in triplicate. The desired in-vial concentrations were  
124 obtained by adding to 1 g ( $\pm 0.0086$  g) of tested matrix 0.1 g ( $\pm 0.0055$  g) of the corresponding standard



125 mixture (approximately 1:10 dilution). The actual weights of matrix and standard mixture were annotated  
126 and used to calculate the in-vial concentrations of the concentration levels.

127 In order to evaluate the matrix effect of the hazelnut paste, the concentration-response curves were  
128 prepared and analyzed in parallel in the two tested matrices: MCT (the model matrix), and the deodorized  
129 roasted hazelnut paste (HZ deo) (matrix-matched blank sample) following the same procedure.

130 All the concentration levels were prepared one day before the HS-GC-IMS measurement, and vortexed for  
131 30 seconds both after preparation and before running the analytical sequence to ensure homogeneity and  
132 analyte mixing in the matrix.

#### 133 **2.4. Identification of alkyl pyrazines in roasted hazelnut paste**

134 Peak identification of alkyl pyrazines in the roasted hazelnut paste samples was performed by comparing the  
135 retention and drift times of pure standard solutions analyzed under the same analytical conditions. The  
136 matching of the experimental drift times (relative to the reactant ion peak - RIP) with those reported in the  
137 built-in library of the VOCal software (see 2.6) was also verified. The identification was confirmed by spiking  
138 increasing concentrations of the target analytes into a hazelnut paste sample.

#### 139 **2.5. Multiple point external standardization**

140 To standardize the IMS detector response and quantify alkyl pyrazines in hazelnut pastes samples an external  
141 calibration approach was applied. A calibration mixture ( $40 \mu\text{g g}^{-1}$  in MCT) was prepared from the stock  
142 solutions of the target analytes identified in the samples (2-methyl-, 2,5-dimethyl-, 2,6-dimethyl-, 2-ethyl-,  
143 and 2,3,5-trimethylpyrazine). The calibration mixture was diluted in the same solvent to obtain the set of  
144 working solutions to prepare the calibration curves. Based on the preliminary results, the deodorized  
145 hazelnut paste was used as a matrix for the external standards to ensure the same matrix effect of the  
146 samples (matrix-matched calibration). The calibration levels were prepared in triplicate directly in the  
147 headspace vials, as described in 2.3, to obtain a set of calibration levels between  $0.01 \mu\text{g g}^{-1}$  and  $1 \mu\text{g g}^{-1}$ .  
148 Blanks were prepared by adding 0.1 g of MCT to 1 g of HZ deo.

#### 149 **2.6. Instrumental set-up**

150 Analyses were carried out on a FlavourSpec GC-IMS system (<sup>3</sup>H-IMS) (G.A.S., Dortmund, Germany) equipped  
151 with a HT2000H headspace autosampler (HTA, Brescia, Italy) and a polar column MXT-Wax 30 m, 0.53 mm  
152 *dc*, 0.5 μm *df* (Restek Corporation, Bellefonte, US). Nitrogen (99.999% purity grade) was used as carrier gas  
153 for the chromatographic separation and as drift gas. The IMS operated in positive ionization mode.  
154 Instrument control and data acquisition were performed with the Sequence Designer software (G.A.S.,  
155 Dortmund, Germany). Data visualization and peak integration were performed with VOCal software (G.A.S.,  
156 Dortmund, Germany).

### 157 2.6.1 HS-GC-IMS parameters

- 158 • GC-IMS

159 Injector and transfer line 1 (inj-oven) temperature 80°C; GC and transfer line 2 (oven-IMS) temperature  
160 60°C, IMS temperature 45°C. GC column flow program: 2 ml/min (constant flow) for 6 min, from min 6  
161 the flow was gradually increased up to 12 ml/min at 16 min, then up to 50 ml/min at 19.5 min, up to 75  
162 ml/min at 22.5 min, up to 124 ml/min at 27 min, up to 150 ml/min at 27 min, ending with 3 min at 150  
163 ml/min constant flow to avoid carryover effects. The total GC runtime was 30 min. IMS drift flow: 150  
164 mL/min.

- 165 • SHS sampling

166 The vials were incubated at 60°C for 20 min under constant agitation. HS syringe temperature was 80°C.  
167 For the first set of experiments (phase 1 - detector response curves and matrix effect assessment) the  
168 injection volume was 0.5 mL. It was increased to 1.0 mL for the quantification (phase 2 - calibration curves  
169 and samples) to improve the detection of alkyl pyrazines present in traces in the samples (2-ethylpyrazine  
170 and 2.6-dimethylpyrazine) or with a higher limit of detection (2,3,5-trimethylpyrazine).

### 171 2.7. Data processing

172 The signal intensity of the selected peaks was extracted as raw data with the quantitative module of the  
173 VOCal software. Peak height was selected as the intensity measure. Microsoft Excel (version 2210, Redmond,  
174 WA, USA) was used to collect data and to perform the calibration and quantification calculations. R.4.2.1. (R

175 Foundation for Statistical Computing, Vienna, Austria) internal statistical functions and the external packages  
 176 “ggplot2” were used to generate the detector response curves.

177


### 178 3) Results and discussion

#### 179 3.1 Evaluation of the GC-IMS response for target alkyl pyrazines

##### 180 3.1.a Detector response curves

181 The IMS detector response for the target alkyl pyrazines was evaluated over 2-3 orders of magnitude of  
 182 concentrations to investigate it both in terms of minimum detectable concentration (low concentration  
 183 levels), and in terms of proton-bound dimer formation (high concentration levels).

184 To identify the suitable concentration range for the detector response curves, a preliminary experiment was  
 185 carried out analyzing solutions of the alkyl pyrazines of interest at 1 and 5  $\mu\text{g g}^{-1}$ . The 2D topographic plot  
 186 (graphical output of the GC-IMS analysis) (**Figure 1**) shows differences in the GC-IMS response for the  
 187 investigated analytes.



188

189 **Fig. 1** - 2D topographic plot (chromatogram x mobilogram) of two solutions of alkyl pyrazines (matrix MCT)  
190 having concentration of  $1 \mu\text{g g}^{-1}$  (left plot) and  $5 \mu\text{g g}^{-1}$  (right plot). 2M: 2-methylpyrazine, 25DM: 2,5-  
191 dimethylpyrazine, 26DM: 2,6-dimethylpyrazine, 2E: 2-ethylpyrazine, TRI: 2,3,5-trimethylpyrazine, TETRA:  
192 2,3,5,6-tetramethylpyrazine, 2E35DM: 2-ethyl-3,5-dimethylpyrazine, 2E36DM: 2-ethyl-3,6-  
193 dimethylpyrazine. M: protonated monomer, D: proton bound dimer.


194 For 2-methyl-, 2,5-dimethyl-, 2,6-dimethyl-, and 2-ethylpyrazine both the protonated monomer and the  
195 proton bound dimer were detectable at  $1 \mu\text{g g}^{-1}$ , and the peak signals were saturated at  $5 \mu\text{g g}^{-1}$ . On the  
196 opposite, for 2,3,5,6-tetramethyl-, 2-ethyl-3,5-dimethyl- and 2-ethyl-3,6-dimethylpyrazine, at  $1 \mu\text{g g}^{-1}$ , only  
197 the protonated monomer peak was observed in the 2D topographic plot, with a low signal intensity. For 2,3,5-  
198 trimethylpyrazine an intermediate behavior was observed. Based on these first observations, the following  
199 concentration ranges were selected to study the detector dynamic range for the target analytes:  $0.1 - 10 \mu\text{g}$   
200  $\text{g}^{-1}$  for 2-methyl-, 2,5-dimethyl-, 2,6-dimethyl-, and 2-ethylpyrazine;  $1 - 100 \mu\text{g g}^{-1}$  for 2,3,5,6-tetramethyl-,  
201 2-ethyl-3,5-dimethyl- and 2-ethyl-3,6-dimethylpyrazine;  $0.1 - 100 \mu\text{g g}^{-1}$  for 2,3,5-trimethylpyrazine.

202  
203 **Figure 2a** reports the concentration response curves of 2-methyl-, 2,5-dimethyl- and 2,3,5-trimethylpyrazine  
204 analyzed in the concentration range  $0.1 - 10 \mu\text{g g}^{-1}$ . Despite their common chemical structure, relevant  
205 differences in the IMS detector response curves can be observed for the alkyl pyrazine homologues. The most  
206 evident difference is the slope of the curves corresponding to the proton bound dimer signals (dark blue),  
207 and, as consequence, the ratio of the two ionized species peak signals. Considering the aromatic ring  
208 substitution, we hypothesized a less energetically favorable formation of the proton bound dimer for higher  
209 homologous, due to the steric hindrance of the methyl groups.

210  
211 The results obtained by comparing the response curve for 2,3,5-trimethylpyrazine within the higher  
212 concentration range ( $1 - 100 \mu\text{g g}^{-1}$ ) with that of 2,3,5,6-tetramethylpyrazine - **Figure 2b** - were consistent  
213 with this hypothesis. While for trimethylpyrazine the signal increase for the proton bound dimer was

214 observable at high concentration, for tetramethylpyrazine, which is more sterically encumbered, the peak  
 215 intensity of the dimer remained low along all the concentration range analyzed.

216 Lastly, as firstly observed in the preliminary experiment, also the lower detectable concentrations, the  
 217 detector sensitivity (curve slope), and the concentration level leading to monomer peak signal saturation  
 218 (curve plateau) were different for the alkyl pyrazines studied.



219

220 **Fig. 2** - SHS-GC-IMS response curves of 2-methylpyrazine (2M), 2,5-dimethylpyrazine (25DM) and 2,3,5-  
 221 trimethylpyrazine (TRI) with concentration range 0.1 – 10  $\mu\text{g g}^{-1}$  (matrix MCT) (A), and 2,3,5-  
 222 trimethylpyrazine (TRI) and 2,3,5,6-tetramethylpyrazine (TETRA) with concentration range 1 – 100  $\mu\text{g g}^{-1}$   
 223 (matrix MCT) (B). Detector response extracted from raw data as peak height [V].


224

225 On the one hand, this experimental evidence proved that the relative intensities of the 2D plot peaks could  
226 not be used as quantitative measure of analyte abundances, because of the remarkable differences in the  
227 GC-IMS response toward different analytes (also for homologous classes having a common chemical  
228 structure, such as alkyl pyrazines). On the other hand, they also highlighted the unsuitability of semi-  
229 quantification based on an internal standard as standardization strategy, because the difference in response  
230 factor of the analytes would lead to inaccurate semi-quantitative results. An external standardization  
231 approach is therefore fundamental to obtain quantitative information for the analytes of interest.

### 232 **3.1.b Matrix effect**

233 Roasted hazelnut paste is a complex matrix in terms of composition, with a high percentage of fat (50-60 %),  
234 but also relevant content of proteins (15%), carbohydrates (17%) and fiber (10%) [23, 24]. An oil-like solvent,  
235 such as the MCT used in this study, could provide a model system for fat-rich food matrices [25], however,  
236 the release of alkyl pyrazines has been reported to be influenced by the protein content in the food matrix  
237 [26].

238 The matrix effect of the hazelnut paste was assessed in the first phase of the research using the deodorized  
239 roasted hazelnut paste to prepare the concentration levels with a matrix matching approach. The comparison  
240 of the IMS detector response curves in the two tested matrices, reported in **Figure 3** (for 2,5-dimethyl- and  
241 2,3,5-trimethylpyrazine), shows a non-negligible matrix effect. In particular, the lower percentage of fat  
242 content of the hazelnut paste, compared to the model matrix MCT, results in higher peak signals for the  
243 target alkyl pyrazines. Based on this result, for the second phase of this research, the deodorized roasted  
244 hazelnut paste (HZ deo) was used to prepare the external standardization calibration levels, in order to  
245 ensure the same matrix effect of the samples.



246

247 **Fig.3** - Comparison of the SHS-GC-IMS response curves (for 2,5-dimethylpyrazine on the left, and 2,3,5-  
 248 trimethylpyrazine on the right) obtained with concentration levels prepared in the model matrix (MCT –  
 249 blue) and with a matrix matching approach (deodorized hazelnut paste - brown). Concentration range 0.1 –  
 250 10  $\mu\text{g g}^{-1}$ . Detector response extracted from raw data as peak height [V].


251

### 252 3.2 Alkyl pyrazines in roasted hazelnut paste samples

253 The following five alkyl pyrazines were identified in the roasted hazelnut paste samples: 2-methyl-, 2,5-  
 254 dimethyl-, 2,6-dimethyl-, 2-ethyl- and 2,3,5-trimethylpyrazine (**Figure 4**). Among them, 2-methylpyrazine and  
 255 2,5-dimethylpyrazine showed the most intense signals in the 2D topographic plot, and both the protonated  
 256 monomer and the proton bound dimer peaks were observed. On the contrary, for 2-ethylpyrazine and 2,3,5-  
 257 trimethylpyrazine only the protonated monomer peak was detected. Lastly, 2,6-dimethylpyrazine monomer  
 258 was observed only for samples corresponding to the longer roasting time (25 min). Under the applied  
 259 analytical conditions, 2,3,5,6-tetramethyl-, 2-ethyl-3,5-dimethyl- and 2-ethyl-3,6-dimethylpyrazine were not  
 260 detected in the roasted hazelnut paste samples.

261





262

263

264 **Fig. 4** - 2D topographic plot of two roasted hazelnut paste samples (Tonda Gentile Romana samples - RO). On  
 265 the left sample obtained with 15 minutes of roasting at 140°C, on the right sample obtained with 25 minutes  
 266 at the same temperature. 2M: 2-methylpyrazine, 25DM: 2,5-dimethylpyrazine, 26DM: 2,6-dimethylpyrazine,  
 267 2E: 2-ethylpyrazine, TRI: 2,3,5-trimethylpyrazine. M: protonated monomer, D: proton bound dimer.

268

269 Based on the results of the first phase of this study, it is clear that the relative intensities of these peaks  
 270 cannot be used to draw quantitative conclusions about the alkyl pyrazine abundances, but a standardization  
 271 approach must be applied.

### 272 3.3 Multiple point external standardization

273 The standardization of the GC-IMS responses for the target alkyl pyrazines was carried out developing a  
 274 quantitative protocol based on multiple point external standardization. Starting from the results of the first  
 275 phase and from literature data, the concentration range selected for the calibration levels was 0.01 - 1 µg  
 276 g<sup>-1</sup>. The ionic species used for the quantification was chosen for each alkyl pyrazine. For 2-methylpyrazine,  
 277 the quantification was carried out on the proton-bound dimer intensity, due to a coelution of the monomer  
 278 with an unknown peak in the samples with higher roasting intensity. In the case of 2,5-dimethylpyrazine



279 calibration both ionic species were used separately. The protonated monomer curve was suitable for the low  
280 concentration range ( $0.01-0.6 \mu\text{g g}^{-1}$ ), while the proton bound dimer curve was applied for the high  
281 concentration range ( $0.2-1 \mu\text{g g}^{-1}$ ). For 2,6-dimethyl-, 2-ethyl- and trimethylpyrazine the protonated  
282 monomer was used, because it was the only detectable species in the hazelnut samples. The linear dynamic  
283 range of GC-IMS is known to be limited to 1-2 order of magnitude [14], and a 2<sup>nd</sup> degree polynomial function  
284 was used to fit the calibration curves of those analytes showing non-linearity of the detector response within  
285 the concentration range (2,6-dimethylpyrazine, 2-ethylpyrazine, and the protonated monomer of 2,5-  
286 dimethylpyrazine). Limit of detection (LODs) and limit of quantification (LOQs) were calculated from the  
287 standard deviation of blanks using the linear approximation of the calibration curves, according to Budzyńska  
288 *et al* [10], and confirmed by analyzing decreasing concentrations of standards. Details of the calibration for  
289 the target alkyl pyrazines are summarized in **Table 1**.

290

	Ion species	Curve fitting	Concentration range ( $\mu\text{g g}^{-1}$ )	Equation	$R^2$	LOD ( $\mu\text{g g}^{-1}$ )	LOQ ( $\mu\text{g g}^{-1}$ )
<b>2M</b>	dimer	linear	0.08-1	$y = 1.418x + 0.0088$	0.9934	0.04	0.08
<b>25DM</b>	monomer	polynomial (2 <sup>nd</sup> )	0.01-0.6	$y = -2.9737x^2 + 3.5177x + 0.1063$	0.9983	< 0.01	0.01
	dimer	linear	0.2-1	$y = 0.647x - 0.0162$	0.9985	0.08	0.2
<b>26DM</b>	monomer	polynomial (2 <sup>nd</sup> )	0.06-1	$y = -0.8385x^2 + 1.4424x + 0.0115$	0.9954	0.03	0.06
<b>2E</b>	monomer	polynomial (2 <sup>nd</sup> )	0.04-1	$y = -0.7912x^2 + 1.4628x + 0.0382$	0.9957	0.02	0.04
<b>TRI</b>	monomer	linear	0.06-1	$y = 0.6693x + 0.0372$	0.9925	0.03	0.06

291

292 **Table 1** - Summary of the multiple point external standardization. 2M: 2-methylpyrazine, 25DM: 2,5-dimethylpyrazine, 26DM: 2,6-dimethylpyrazine, 2E: 2-  
 293 ethylpyrazine, TRI: 2,3,5-trimethylpyrazine.  $R^2$ : regression coefficient; LOD: limit of detection; LOQ: limit of quantitation; reported in  $\mu\text{g g}^{-1}$ .

294 **3.4 Determination of alkyl pyrazines in roasted hazelnut paste samples**

295 The external calibration developed to standardize the GC-IMS response was applied to determine the content  
296 of the identified alkyl pyrazines in the roasted hazelnut paste samples.

297 As mentioned in paragraph 3.2, 2,6-dimethylpyrazine monomer was detectable only for samples  
298 corresponding to higher roasting intensity (roasting time 25 min). Even in these analyses, however, the peak  
299 intensity was lower than the LOQ. Therefore, it was not possible to determine the content of 2,6-  
300 dimethylpyrazine in the samples.

301 The results obtained for other analytes are reported in **Table 2**. As expected, an increase in the alkyl pyrazines  
302 content according to the roasting times was observed. At the same time, TGL samples showed an overall  
303 higher amount of all the pyrazines quantified, while lower amounts were observed for AK samples. 2-  
304 methylpyrazine and 2,5-dimethylpyrazine are the most abundant alkyl pyrazines in all the samples analyzed,  
305 while 2-ethyl pyrazine is the less abundant.

306 2,3,5-trimethyl pyrazine is detectable in the 2-D plot as low intensity monomer peak, therefore observations  
307 based on the raw signal intensities would have led to a considerable underestimation of the amount of this  
308 analyte compared to 2-methyl and 2,5-dimethylpyrazine (both monomer and dimer detectable). On the  
309 contrary, the results obtained applying the response standardization show that the actual amount of 2,3,5-  
310 trimethylpyrazine is comparable to those of 2-methylpyrazine in many samples (low roasting intensity  
311 samples), despite the differences in the signal patterns. The standardization enables, therefore, comparative  
312 evaluation among analytes within the same sample.

313 Lastly, after the standardization, a single analytical output is obtained (the content expressed in  $\mu\text{g g}^{-1}$ ),  
314 instead of multiple peak intensity values for the same analyte, making the results more easily interpretable.

	TGL			RO			AK		
	15'	20'	25'	15'	20'	25'	15'	20'	25'
<b>2M</b>	0.18 ± 0.01	0.36 ± 0.00	0.59 ± 0.01	0.08 ± 0.01	0.28 ± 0.00	0.55 ± 0.01	< LOQ	0.19 ± 0.00	0.42 ± 0.01
<b>25DM</b>	<i>0.30 ± 0.02</i>	0.49 ± 0.00	0.65 ± 0.01	<i>0.26 ± 0.00</i>	0.36 ± 0.02	0.53 ± 0.00	<i>0.12 ± 0.00</i>	<i>0.24 ± 0.00</i>	0.40 ± 0.00
<b>2E</b>	< LOQ	0.07 ± 0.01	0.11 ± 0.01	< LOQ	0.06 ± 0.01	0.09 ± 0.01	< LOQ	0.05 ± 0.00	0.09 ± 0.01
<b>TRI</b>	0.17 ± 0.01	0.27 ± 0.01	0.35 ± 0.01	0.09 ± 0.01	0.21 ± 0.01	0.30 ± 0.01	0.07 ± 0.00	0.14 ± 0.01	0.24 ± 0.00

315

316 **Table 2** - Quantitative results of the analysis of alkyl pyrazines in roasted hazelnut paste samples (different roasting times: 15 – 20 – 25 minutes) by SHS-GC-IMS.317 Quantitative values, expressed as the mean of three replicates ± standard deviation, reported in  $\mu\text{g g}^{-1}$ . For 25DM values obtained from protonated monomer

318 calibration are reported in italics. 2M: 2-methylpyrazine, 25DM: 2,5-dimethylpyrazine, 26DM: 2,6-dimethylpyrazine, 2E: 2-ethylpyrazine, TRI: 2,3,5-

319 trimethylpyrazine. **TGL**: Tonda Gentile delle Langhe, **RO**: Tonda Gentile Romana, **AK**: Akçakoca

#### 4) Conclusions

320

321 SHS-GC-IMS is becoming widely applied for the analysis of food volatile organic compounds. In this study we  
322 evaluated its application within the context of industrial process monitoring, in particular to determine the  
323 content of alkyl pyrazines in roasted hazelnut pastes. When IMS is operated in positive ionization mode it  
324 works as a non-selective detector and the volatile compounds with proton affinity higher than water are  
325 ionized and detected. In the first phase of the study, relevant differences in IMS response were observed for  
326 alkyl pyrazine homologous both in terms of minimum detectable concentration and sensitivity, but also  
327 concerning the proton-bound dimer formation. These differences highlighted that a response  
328 standardization approach is fundamental for IMS detection when comparative evaluations among analyte  
329 abundances are performed.

330 The approach proposed in the second section of this article is based on multiple point external  
331 standardization and was carried out applying in-matrix calibration to avoid the matrix effect, which was  
332 observed for the hazelnut paste. It allowed to obtain a single standardized analytical output and to compare  
333 the abundances of alkyl pyrazines in hazelnut paste samples from different geographical origins and  
334 processed with different roasting intensity.

335 In conclusion, HS-GC-IMS is suitable as a rapid analytical platform to monitor VOCs during industrial  
336 processing, however, before any quantitative assumptions are drawn, the IMS response needs to be carefully  
337 evaluated and standardized for the analytes of interest, even in the case of homologous compounds.

338

339 **References**

- 340 [1] G. A. Eiceman, Z. Karpas, and H. H. J. Hill, *Ion Mobility Spectrometry*, 3rd Edition. Taylor&Francis  
341 Group, 2013.
- 342 [2] G. A. Eiceman, "Advances in Ion Mobility Spectrometry: 1980—1990," *Crit. Rev. Anal. Chem.*, vol. 22,  
343 no. 1–2, pp. 471–490, 1991, doi: 10.1080/10408349108055024.
- 344 [3] R. H. Louis, H. H. Hill, and G. A. Eiceman, "Ion Mobility Spectrometry in Analytical Chemistry," *Crit.*  
345 *Rev. Anal. Chem.*, vol. 21, no. 5, pp. 321–355, 1990, doi: 10.1080/10408349008050848.
- 346 [4] M. J. Cohen and F. W. Karasek, "Plasma chromatography - a new dimension for gas chromatography  
347 and mass spectrometry," *J. Chromatogr. Sci.*, vol. 8, no. 6, 1970, doi: 10.1093/chromsci/8.6.330.
- 348 [5] S. Wang, H. Chen, and B. Sun, "Recent progress in food flavor analysis using gas chromatography–  
349 ion mobility spectrometry (GC–IMS)," *Food Chem.*, vol. 315, no. January, p. 126158, 2020, doi:  
350 10.1016/j.foodchem.2019.126158.
- 351 [6] S. Gu, J. Zhang, J. Wang, X. Wang, and D. Du, "Recent development of HS-GC-IMS technology in rapid  
352 and non-destructive detection of quality and contamination in agri-food products," *TrAC - Trends in*  
353 *Analytical Chemistry*, vol. 144. 2021, doi: 10.1016/j.trac.2021.116435.
- 354 [7] D. Cavanna, S. Zanardi, C. Dall'Asta, and M. Suman, "Ion mobility spectrometry coupled to gas  
355 chromatography: A rapid tool to assess eggs freshness," *Food Chem.*, vol. 271, no. July 2018, pp.  
356 691–696, 2019, doi: 10.1016/j.foodchem.2018.07.204.
- 357 [8] N. Jurado-Campos, A. Martín-Gómez, D. Saavedra, and L. Arce, "Usage considerations for  
358 headspace-gas chromatography-ion mobility spectrometry as a suitable technique for qualitative  
359 analysis in a routine lab," *J. Chromatogr. A*, vol. 1640, p. 461937, 2021, doi:  
360 10.1016/j.chroma.2021.461937.
- 361 [9] W. Zhu, F. Benkwitz, B. Sarmadi, and P. A. Kilmartin, "Validation Study on the Simultaneous  
362 Quantitation of Multiple Wine Aroma Compounds with Static Headspace-Gas Chromatography-Ion

- 363 Mobility Spectrometry," *J. Agric. Food Chem.*, 2021, doi: 10.1021/acs.jafc.1c06411.
- 364 [10] E. Budzyńska, S. Sielemann, J. Puton, and A. L. R. M. Surminski, "Analysis of e-liquids for electronic  
365 cigarettes using GC-IMS/MS with headspace sampling," *Talanta*, vol. 209, no. September 2019, p.  
366 120594, 2020, doi: 10.1016/j.talanta.2019.120594.
- 367 [11] C. J. Denawaka, I. A. Fowles, and J. R. Dean, "Evaluation and application of static headspace-  
368 multicapillary column-gas chromatography-ion mobility spectrometry for complex sample analysis,"  
369 *J. Chromatogr. A*, vol. 1338, 2014, doi: 10.1016/j.chroma.2014.02.047.
- 370 [12] R. Brendel, S. Schwolow, S. Rohn, and P. Weller, "Comparison of PLSR, MCR-ALS and Kernel-PLSR for  
371 the quantification of allergenic fragrance compounds in complex cosmetic products based on  
372 nonlinear 2D GC-IMS data," *Chemom. Intell. Lab. Syst.*, vol. 205, no. May, p. 104128, 2020, doi:  
373 10.1016/j.chemolab.2020.104128.
- 374 [13] "https://webbook.nist.gov/cgi/cbook.cgi?ID=C7732185&Mask=20."  
375 https://webbook.nist.gov/cgi/cbook.cgi?ID=C7732185&Mask=20 (accessed Feb. 21, 2023).
- 376 [14] M. J. Cardador, N. Jurado-Campos, and L. Arce, "Ion mobility detectors for gas chromatography," in  
377 *Gas Chromatography*, 2021.
- 378 [15] N. Jurado-Campos, R. Garrido-Delgado, B. Martínez-Haya, G. A. Eiceman, and L. Arce, "Stability of  
379 proton-bound clusters of alkyl alcohols, aldehydes and ketones in Ion Mobility Spectrometry,"  
380 *Talanta*, vol. 185, no. March, pp. 299–308, 2018, doi: 10.1016/j.talanta.2018.03.030.
- 381 [16] C. Wang, Z. Zhang, X. Zhang, X. Tian, K. Chen, and X. Zeng, "Characterization of Volatile Compounds  
382 by HS-GC-IMS and Chemical Composition Analysis of Colored Highland Barley Roasted at Different  
383 Temperatures," 2022.
- 384 [17] S. Ge, Y. Chen, S. Ding, H. Zhou, and L. Jiang, "Changes in volatile flavor compounds of peppers  
385 during hot air drying process based on headspace-gas chromatography-ion mobility spectrometry (HS-GC-IMS)," no. December 2019, 2020, doi: 10.1002/jsfa.10341.
- 386

- 387 [18] R. Müller and S. Rappert, "Pyrazines: Occurrence, formation and biodegradation," *Applied*  
388 *Microbiology and Biotechnology*, vol. 85, no. 5. 2010, doi: 10.1007/s00253-009-2362-4.
- 389 [19] N. M. Fayek, J. Xiao, and M. A. Farag, "A multifunctional study of naturally occurring pyrazines in  
390 biological systems; formation mechanisms, metabolism, food applications and functional  
391 properties," *Critical Reviews in Food Science and Nutrition*. 2021, doi:  
392 10.1080/10408398.2021.2017260.
- 393 [20] J. García-Lomillo and M. L. González-SanJosé, "Pyrazines in thermally treated foods," in *Encyclopedia*  
394 *of Food Chemistry*, 2018.
- 395 [21] L. Hashim and H. Chaveron, "Use of methylpyrazine ratios to monitor the coffee roasting," *Food Res.*  
396 *Int.*, vol. 28, no. 6, pp. 619–623, 1995, doi: 10.1016/0963-9969(95)00037-2.
- 397 [22] L. Nicolotti, C. Cordero, C. Bicchi, P. Rubiolo, B. Sgorbini, and E. Liberto, "Volatile profiling of high  
398 quality hazelnuts (*Corylus avellana* L.): Chemical indices of roasting," *Food Chem.*, vol. 138, no. 2–3,  
399 pp. 1723–1733, 2013, doi: 10.1016/j.foodchem.2012.11.086.
- 400 [23] "[https://fdc.nal.usda.gov/fdc-app.html#/food-details/2342979/nutrients.](https://fdc.nal.usda.gov/fdc-app.html#/food-details/2342979/nutrients)"  
401 <https://fdc.nal.usda.gov/fdc-app.html#/food-details/2342979/nutrients> (accessed Feb. 21, 2023).
- 402 [24] K. Król and M. Gantner, "Morphological traits and chemical composition of hazelnut from different  
403 geographical," *Agric.*, vol. 10, no. 9, pp. 1–15, 2020, doi: 10.3390/agriculture10090375.
- 404 [25] R. G. Buttery, D. G. Guadagni, and L. C. Ling, "Flavor Compounds: Volatilities in Vegetable Oil and Oil-  
405 Water Mixtures. Estimation of Odor Thresholds," *J. Agric. Food Chem.*, vol. 21, no. 2, 1973, doi:  
406 10.1021/jf60186a029.
- 407 [26] Y. J. Ma *et al.*, "Effect of alkyl distribution in pyrazine on pyrazine flavor release in bovine serum  
408 albumin solution," *RSC Adv.*, vol. 9, no. 63, 2019, doi: 10.1039/c9ra06720e.

409



## Highlights

- The determination of alkyl pyrazines in roasted hazelnut pastes by GC-IMS is presented.
- The SHS-GC-IMS response curves of alkyl pyrazine homologous differ significantly.
- Detectable concentrations, sensitivity, proton-bound dimer formation are different.
- Abundance comparison among analytes requires a response standardization.
- A matrix-matched multiple point external standardization was applied.

**Declaration of interests**

The authors declare that they have no known competing financial interests or personal relationships that could have appeared to influence the work reported in this paper.

The authors declare the following financial interests/personal relationships which may be considered as potential competing interests:

Journal Pre-proof